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Read the operating instructions carefully before using the appliance for the first time and preserve this booklet for later reference. Pass this booklet on to whoever might acquire the appliance at a later date.

Danger of electric shock!

- Connect the appliance only to correctly installed and properly earthed power sockets. Ensure that the rating of the local power supply tallies completely with the details given on the rating plate of the appliance.
- Always remove the plug from the wall socket when you attach accessories, in the case of malfunctions, before cleaning the appliance and when it is no longer required for use.
- NEVER submerge the appliance in liquids.

Should liquids nonetheless permeate into the appliance housing, immediately remove the plug from the wall socket. Arrange for it to be repaired by a qualified specialist.

DO NOT use the appliance in the bath, in the vicinity of water-filled baths, swimming pools or washbasins.

- NEVER grasp the appliance, power cable or power plug with wet hands.
 - Always remove the power cable from the socket by pulling on the plug, do not pull on the cable itself.
- Do not hold or carry the appliance with the power cable.

- Always remove the power cable from the socket by pulling on the plug, do not pull on the cable itself. Do not hold or carry the appliance with the power cable.

Risk of Fire!

- Do not leave the appliance unattended when in use.
- To avoid a build-up of heat, NEVER cover the appliance or the ventilation openings when it is in use. Do not use the appliance under a blanket or at other poorly ventilated locations.
- Avoid the addition of heat from other sources e.g. from direct sunlight, heating appliances, other devices etc.!
- Do not use an extension cable. In the event of danger, the plug must be quickly and easily accessible.

Risk of personal injuries!

- This appliance is not suitable for children or for people whose physical, sensorial or intellectual faculties prevent the safe usage of the appliance without support or supervision. Children should be watched-over to ensure that they do not play with the appliance.
- Lay the power cable in such a way that no one can tread on or trip over it.

Medical Warning Note

This Massage Appliance is intended for personal use and is not a substitute for physiotherapeutic or medicinal treatments.

The appliance may not be used ...

- by people fitted with heart pace-makers
- by people with heart problems or pains of unknown origin
- in the event of feverish colds, varicose veins, thrombosis, phlebitis, jaundice, diabetes, nerve ailments such as sciatica, acute inflammations or during pregnancy.
- on swollen, inflamed or injured body parts.
- for the massage of bones, especially the spine and the skull.
- close to the eyes or other sensitive body parts.
- on small children.

If you are unsure about your current medical condition, obtain advice from your doctor before using the Massage Appliance.

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